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Produktname

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Instructions for use

Product name

Stethoscope Rappaport, with accessory

Purpose

The stethoscope is a medical diagnostic instrument for auscultating various body sounds.

Product description

The stethoscope is a diagnostic instrument and can be used to auscultate heart, lung, intestinal and other body sounds. Depending on the examination, the stethoscope is to be placed on the desired body region and used under slight pressure.

Properties:

The stethoscope is adjustable in 5 different ways. There are 2 diaphragms (small & large) and 3 chest pieces in different sizes supplied as a set. Each specific Purpose thus receives the corresponding membrane or bell.

1. **Large diaphragm:** This chest piece is used to capture low frequency and diastolic heart sounds. In addition, the chest piece is used to detect first, second and third heart sounds. It is also used to detect high-frequency heart sounds.
2. **Small diaphragm:** Used to detect high-frequency heart sounds in children.
3. **Bell size for adults:** This bell has its special suitability in the determination of low and medium frequency tones and sounds.
4. **Medium bell size:** This bell is used to perceive low and medium frequency sounds and noises, especially in narrow places, such as between the ribs.
5. **Bell size for children:** This bell is specifically for the examination of infants and small children. Sounds and sounds of children's hearts, in the low and medium frequency range, are best recorded best with this bell size.

Handling:

1. Place the earpieces in your ears.
2. Test whether the diaphragm is in the working position by tapping it lightly with your finger.
3. If no knocking sound is audible, check that the stethoscope is correctly assembled.
4. Tap again to check that you hear sounds, and that the stethoscope is now correctly set for use.
5. Depending on the examination, the stethoscope is to be placed on the desired body region and used under slight pressure.
6. When using a different chest piece, please proceed as described under point 1-4.

Cleaning and disinfection:

- Do not use oils and harsh cleaning agents.
- The entire stethoscope can be wiped with a damp cloth.
- To disinfect, wipe the stethoscope with a suitable disinfectant.
- To clean the earpieces, remove them from the ear hook.
- The stethoscope must not be immersed in liquids or autoclaved.
- Make sure that the stethoscope is dry before assembling it.

Durability

The stethoscope must not be used beyond the expiry date or if premature signs of material wear are visible.

Storage instructions

Store in a dry place between 10 - 35°C. Protect from direct sunlight.

Safety instructions:

- Do not use if the packaging or the product is damaged.
- Disinfect the chestpiece before and after use.
- After use, put the stethoscope back in its packaging to avoid damaging the chestpiece.
- Avoid putting pressure on the chestpiece as it could be damaged.
- If the chestpiece or tube is damaged, the stethoscope must no longer be used. Errors in auscultation may occur.
- Avoid hitting the chestpiece when the earpieces are in the ear to protect the eardrum.
- The stethoscope must not be immersed in liquids.
- If there is limited functionality after cleaning, liquid may have penetrated, or the device may have been improperly assembled. After sufficient drying time, reassemble the device and carry out the test as described under Handling. If the device does not function after this, it is damaged and must be disposed of.
- During drying, the stethoscope must not be exposed to increased temperatures. Drying should be done at room temperature.

Legal note

Please report serious incidents that have occurred in connection with this product to the manufacturer and the competent authority of the Member State in which the user and/or patient is located.

Disposal instructions

The stethoscope must be disposed of in accordance with the applicable regulations.

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BODY PRODUCTS RELAX Pharma und Kosmetik GmbH