

SPECIFICATIONS

ROMED PILL PULVERIZER

REF: PP-400



Medical device conform MDR (EU) 2017/745

Pill pulverizer can be used to pulverize pills for easier swallowing or to ingest part of the dose of the medication. It can also serve as a small storage for medicines in the blue part and upper part can be used as a pill cutter. It is safe and easy to use.

With locking button, protect user from potential injury from sharp stainless steel blade.

Materials: Plastic, stainless steel

Height: 6 cm

Diameter: 4 cm

Weight: 46 g

DIRECTIONS FOR USE

1. To cut pills, simple open the white cap.
2. Place pill into the 'V' shaped holder.
3. Close cover to split pill.
4. Repeat opening process and remove pill halves.
5. Lift the white part and store the pill halves in the blue part with the separated sides
6. To store medication, lift the blue part with the separated sides
7. To pulverize pills unscrew the blue part from the white bottom part, put in pill and screw blue part on white part, then unscrew again and pill will be pulverized
8. Close cover and store in closed position.

The pill pulverizes do not contain latex or phthalates.

Packing: per piece in a box, 20 pieces in an inner box and 20 inner boxes in an export carton.