FIRST AID MEASURES

CONTENT LIST DIN 13157 & DIN 13164-2022

NL: Stappenplan Eerste Hulp | Inhoud
DE: Erste Hilfe Maßnahmen | Inhalt
FR: Mesures de premiers secours | Contenu
ES: Medidas de primeros auxilios | Contenido
PT: Medidas de primeiros socorros | Conteúdo
IT: Misure di primo soccorso | Contenuto
HU: Elsősegélynyújtási intézkedések | Tartalom
RS: Mere prve pomoći | Садржај
PL: Środki pierwszej pomocy | Zawartość
RO: Măsurile de prim ajutor | Conținut
GR: Μέτρα πρώτων βοηθειών | Περιεχόμενα
BS: Mjere prve pomoći | Sadržaj
UA: Заходи першої допомоги | Вміст

ROMED - HOLLAND
VAN OOSTVEEN MEDICAL B.V.
HERENWEG 269, 3648 CH WILNIS, THE NETHERLANDS
WWW.ROMED.NL

CALL 1-1-2 	✓ Ensure a safe environment ✓ Inspect the injured person ✓ Call for help ✓ Provide first aid and care as best you can	TRANSPORT ONLY IN EMERGENCY SITUATIONS: If it is absolutely necessary to move the injured person for safety reasons.			
	 Carefully lift the person's shoulders	 lift the head and	 Place the injured person in a sitting position	 Put your hands under their armpits and firmly grasp one of their forearms	 Lift the injured person and move them backwards.

				CALL FOR HELP Where? Give exact details of the location of the emergency. What? What kind of help do you need? What happened? Give a brief description of the accident. How many? Indicate how many people are injured. What kind of injuries? Indicate whether there are any life-threatening injuries. Wait! The operator will end the call once they have all the necessary information.
Switch on the car's hazard warning lights	Ensure you have a powder extinguisher	Always wear a safety vest	Place a warning triangle at a safe distance Motorway 100 m Ordinary road 50m	

CHECK THE INJURED PERSON:

1. Find out what happened 2. Reassure the victim 3. Help the victim as much as possible where they are lying

Check consciousness ✓ Talk to the victim and gently shake their shoulders		Check breathing ✓ Check if the chest is moving ✓ Feel for airflow with your cheek ✓ Listen for breathing sounds		30 Chest compressions – 2 Breaths	
IF NOT BREATHING, CALL 1-1-2, Start CPR					
1. Push down in the middle of the chest, 2. Push down about 5-6 cm, about 100-120 times per minute, 3. Take a normal breath and place your lips around the person's mouth, 4. Blow until you see the chest rise 5. Continue this procedure until the AED is ready and indicates that you should stop temporarily.					
METHOD OF ARTIFICIAL RESPIRATION MOUTH TO MOUTH					
Lay the unconscious victim on their back, with their head tilted back as far as possible	Kneel next to the unconscious person's head, with the palm of your hand on the tip of their chin; push the lower jaw forward and upward against the upper jaw.	With your other free hand, pinch the victim's nose closed and blow air into their mouth. Adults: 12 breaths per minute. Children: 20 breaths per minute	This method can also be used without direct mouth-to-mouth contact, air can be blown into the victim's mouth using a mouthpiece. In this case, the fingers of the first hand that are still free (see phase 2) must be used to firmly press the victim's lips closed around the mouthpiece used. fingers of the first hand (see phase 2) must firmly press the victim's lips closed around the mouthpiece used		

MINISTERIAL MEMORANDUM AND INSTRUCTIONS FOR USE OF THE FIRST AID KIT

First aid before the arrival of the doctor

I. NON-BLEEDING INJURIES

Wash with soapy water and apply a sterile gauze pad with a bandage to the wound.

II. SKIN ABRASIONS OR MINOR INJURIES

If there is an abrasion or a small opening in the epidermis with insignificant blood loss, after washing and drying the injured area, apply a small cotton wool pad soaked in iodine spirit to the injury. cover the injured area with a sterile gauze dressing and a bandage. If the abrasion or small wound is on the fingers or toes, instead of using a dressing with iodine-soaked cotton wool, wrap the wound once or twice with a sterile adhesive dressing.

Abrasion: Clean the abrasion with tap water, using a soft brush if necessary. Remove any remaining stones with sterile tweezers and/or sterile gauze. Disinfect the graze with a bactericidal agent. Covering the graze is only necessary if the clothing is abrasive.

III. BLEEDING INJURIES (General instructions)

It is forbidden to wash the wound, touch it with your fingers or bring it into contact with any object other than a sterile dressing.

A. Simple wounds or bruises (without excessive blood loss): Unroll the bandage without touching the gauze dressing. Place the gauze dressing on the wound and secure it with the bandage.

Simple wound* (little blood): Remove all dirt from the wound with running lukewarm tap water. Allow to bleed for a while. Pat the skin around the wound dry with a sterile gauze pad. Then disinfect the skin around the wound with a bactericidal agent. Cover (with a plaster or other dressing)

B. Wounds with excessive bleeding:

Serious wound*: Reduce bleeding by keeping the injured body part higher than the heart. Cover the wound with a sterile quick dressing. Consult a doctor.

Serious external bleeding*: Reduce bleeding by raising the injured body part and applying pressure to the wound and/or the supplying artery. Apply a pressure bandage (e.g. a quick bandage wrapped in synthetic cotton wool and a sturdy bandage or cloth to apply pressure to the wound). Rest and support the injured body part. Go to a doctor.

1. Bleeding without a jet: Unroll a bandage and press it onto the wound. If necessary, place two or three on top of each other.

2. Bleeding from a limb artery or vein (Intermittent spurting of bright red blood from the wound = arterial wound. Loss of black blood = vein wound).

a) Never remove blood clots that have formed in the wound or on clothing. b) Keep the injured limb raised.

c) Apply pressure to the wound using a loose but unrolled bandage and apply it to the wound. d) Secure this pad with an unrolled bandage. e) If bleeding continues, add another pad and continue to apply pressure to the wound. f) If necessary, tie off the artery.

IV. BRUISES, SPRAINS, BROKEN BONES, DISLOCATIONS

Cool with water, ice or a cold pack for at least 10 minutes. Never allow water to come directly onto the wound and always place a cloth between the skin and the ice or cold pack.

ICE rule:

I Immobilise. Ensure that the affected body part does not move or is not used to support itself.

C Compression. Apply pressure to the affected body part to prevent further swelling. Preferably apply a pressure bandage consisting of synthetic cotton wool and an elastic bandage or cloth.

E Elevation. Raise the affected body part (leg up or arm in a sling). Be aware of the possibility of a broken bone.

Broken bone: Keep as immobile as possible. Provide support and rest. Cover open fractures with a sterile dressing. Call for expert help.

Do not attempt to straighten broken bones or dislocated limbs. Do not remove clothing. Bandage the broken limb and immobilise it. In case of a wound: Do not tear off clothing. Cut or tear them away so that the wound is exposed. Apply one or more bandages as indicated in III. Immobilise the limb.

V. BURNS (General provisions)

Burn: Cool immediately (for at least 10 minutes). Preferably with running lukewarm water from the tap.

Apply a sterile dressing to second-degree (red with blisters) or third-degree (grey-white, boiled or blackened) burns using a metaline dressing. Rest and support the injured body part. Consult a doctor for extensive first-degree, second-degree and third-degree burns. Anyone who has suffered a burn should be treated as an injured person, i.e. they should receive first aid at the scene and then be taken to the first-aid post as quickly as possible and under the best conditions.

A. Minor burns: Apply a bandage as indicated in III. § A.

B. Severe burns: 1. Caused by fire, Never tear off clothing. Cover the exposed areas with bandages. Wrap the victim in a warm bed sheet. Warm the victim with blankets and keep them out of the cold. 2. By acids, Remove clothing soaked with corrosive substances. Sprinkle the affected areas with aseptic chalk powder.

3. Alkali (potassium, soda, etc.) Remove clothing soaked with corrosive substances. Apply compresses moistened with water mixed with vinegar (one bottle of acetic acid to one litre of boiled water). 4. For burns to the face and eyes caused by fire or chemicals, refrain from touching the wounds in any way. These wounds require careful dressing and immediate transfer of the affected person to a place that is fully equipped for appropriate treatment.

VI. FAINTING, CHOKING, DROWNING, ELECTROCUTION

Fainting: Lay the victim down and allow them to remain lying down for at least 10 minutes after regaining consciousness. In the meantime, monitor breathing. If you are unsure of the cause, consult a doctor.

1. Lay the affected person flat on their back in the open air.

2. Then first make sure that the openings of the respiratory organs, namely the nostrils, mouth and throat, are completely clear.

To this end: a) Blow gently into the nostrils of the affected person to clear them of any mucus or dust that may have accumulated there; b) Open the mouth and hold the jawbones apart with a cork or piece of wood; c) Grasp the tongue with a cloth (handkerchief or other) and pull it out, in order to remove mucus and other substances from the mouth and induce vomiting; d) Quickly clean the mouth, nose and throat with your finger. It is advisable to hold drowning victims and those who have been suffocated by dust with their head downwards for a few moments, as the mucus and dust collected in the lung tubes can be removed more easily in this unusual position.

3. Facilitate breathing: To this end: a) Clear the neck of any obstructions (collar, tie, scarf, etc.);

b) Open the clothing wide (jacket, waistcoat, shirt, trousers) and remove any belts or straps in order to free the chest and abdomen of all obstructions and facilitate the action of the respiratory muscles.

4. If the three aforementioned measures are ineffective, artificial respiration should be administered as soon as possible.

5. Continue artificial respiration until breathing movements reappear.

6. Do not give up, continue for hours if necessary and in any case until the arrival of a doctor.

7. If help is available, rub the body and limbs of the affected person vigorously and slap their face with a slightly damp cloth.

8. When the affected person begins to breathe, continue to provide care. They will still require special assistance: a) Warm the affected person with blankets and vigorous rubbing;

b) Pay attention to breathing so that artificial respiration can be applied immediately if spontaneous breathing stops.

VII. SUDDEN HEAT

Remove clothing. Wash the entire body thoroughly with cold water. If necessary, apply artificial respiration.

VIII. POISONING

1. Induce vomiting by tickling the throat with a finger, a feather, a brush, etc. Resume tickling after giving the patient plenty of salt water to drink. In addition.

2. In case of poisoning by an acid: Administer chalk powder stirred in water (one tablespoon of chalk powder for every 250 grams of water). 3. In case of poisoning by an alkali (potassium, soda, etc.): Administer water mixed with vinegar (6 to 8 tablespoons of vinegar per litre of water). 4. In all cases of poisoning, seek medical assistance immediately and, while waiting, prepare the following items, which are needed for gastric lavage: a syringe, a funnel, plenty of lukewarm water and a little olive oil.

NOTE: *Tetanus can be contracted from a non-surgical wound that has occurred under non-sterile conditions. The tetanus bacterium is found mainly in soil, street dirt and the faeces of herbivorous animals. Consult a doctor within 48 hours if your last tetanus vaccination was more than 10 years ago or if you are unsure

INSTRUCTIONS

QUICK BANDAGES: These consist of a cotton wool or gauze pad to be placed on the wound and a band to hold the pad in place. If the band is tightened sufficiently, a pressure bandage can easily be applied to stop bleeding.

CAMBIC OR GAUZE BANDAGE: To secure a bandage or to dress a wound that is almost healed.

CUTTING BANDAGE: To be used for sprains and dislocations or to apply a pressure and elastic bandage. In other words, to support a limb without immobilising it immobile. Adhesive rubber plaster with zinc oxide to secure the bandages. If it sticks too firmly, remove with a little ether.

TOXIC BLOODSTOPPER: To stop significant blood flow. Place between the wounded vein and the heart. Do not leave the tourniquet on for longer than 20 minutes. If necessary, loosen it slightly and then reapply. Indicate the exact time of application for the doctor's information.

IODINE ALCOHOL: Disinfectant for bruises and minor wounds.

ETHER: For degreasing and cleaning the skin and wounds.

FIRST AID KIT DIN 13164	FIRST AID KIT DIN 13157
Specification contents	Specification contents
1 adhesive tape 2.5cmx500cm (DIN 13019 - A5x2.5)	1 adhesive tape 2.5cmx500cm (DIN 13019 - A5x2.5)

FIRST AID KIT DIN 13164	FIRST AID KIT DIN 13157
4 adhesive plasters (DIN 13019 -E10x6) 2 fingertip plasters (DIN13019-5x4) 2 fingertip plasters (DIN 13019 -12x2) 2 adhesive plasters (DIN 13019 - 7,2x1,9) 4 adhesive plasters (DIN 13019 - 7,2x2,S) 1 dressing pack 6cmx8cm STERILE (DIN 13151-K) 2 dressing pack 8cmx10cm STERILE (DIN 13151-M) 1 dressing pack 10cmx12cm STERILE (DIN 13151-G) 1 surgical sheet 60cm x 80cm STERILE (DIN 13152-A) 3 elastic bandages 8cm x4m (DIN 61631-MB-BCV /CO) 2 elastic bandages 6cm x4m (DIN 61631-MB-6CV /CO) 1 rescue cover 210cm x 160 cm 6 wound dressing 10cm x 10cm STERILE 1 triangularbandage (DIN 13168-D) 1 pair of scissors (DIN 58279 -A145) 4 gloves in polybag (DIN EN 455-1 & DIN EN 455-2) 2 cleansing wipes 2 face masks (DIN EN 14683) 1 first Aid instruction booklet 1 list of contents	12 adhesive plasters (DIN 13019 - E10x6) 6 fingertip plasters (DIN13019-Sx4) 6 fingertip plasters (DIN 13019 -12x2) 6 adhesive plasters (DIN 13019 - 7,2x1,9) 12 adhesive plasters (DIN 13019 - 7,2x2,5) 1 dressing pack 6cmx8cm STERILE (DIN 13151-K) 3 dressing pack 8cmx10cm STERILE (DIN 13151-M) 1 dressing pack 10cmx12cm STERILE (DIN 13151-G) 1 surgical sheet 60cm x 80cm STERILE (DIN 13152-A) 6 wound dressing 10cm x 10cm STERILE 2 eye pads STERILE 1 instant cold packs 1 rescue cover 210cm x 160 cm 2 elastic bandages 8cm x 4m (DIN 61631-MB-8CV /CO) 2 elastic bandages 6cm x 4m (DIN 61631-MB-6CV /CO) 2 triangular bandages (DIN 13168-D) 1 pair of scissors (DIN 58279 -8190) 2 PE-bags 30cm x 40cm 5 blanket made of non-woven fabric 2 cleansing wipes 4 gloves in polybag (DIN EN 455-1 & DIN EN 455-2) 2 face masks (DIN EN 14683) 1 first aid instructions 1 list of contents

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