ZINKGELBINDE / NOBAZINCULAST® / NOBAVARO®

- D** **Gebrauchsinformation**
ZINKGELBINDE / NOBAZINCULAST® / NOBAVARO®
- GB** **Information for Use**
ZINKGELBINDE / NOBAZINCULAST® / NOBAVARO®
- NL** **Gebruiksaanwijzing**
ZINKGELVERBAND / NOBAZINCULAST® / NOBAVARO®
- DK** **Brugsinformation**
ZINKGELFORBINDING / NOBAZINCULAST® / NOBAVARO®
- F** **Notice d'utilisation**
BANDE AU GEL DE ZINC / NOBAZINCULAST® / NOBAVARO®
- I** **Informazioni per l'uso**
BENDA AL GEL DI ZINCO / NOBAZINCULAST® / NOBAVARO®
- PL** **Informacje dotyczące użytkowania**
ŻELOWEGO BANDAŻA CYNKOWEGO / NOBAZINCULAST® / NOBAVARO®
- ES** **Información de uso**
VENDA DE GEL DE ZINC / NOBAZINCULAST® / NOBAVARO®
- PT** **Informações de uso**
LIGADURA DE GEL DE ZINCO / NOBAZINCULAST® / NOBAVARO®
- SV** **Bruksanvisning**
ZINKGELFÖRBAND / NOBAZINCULAST® / NOBAVARO®
- CZ** **Informace k použití**
GELOVÉ OBVAZY SE ZINKEM / NOBAZINCULAST® / NOBAVARO®
- LT** **Naudojimo informacija**
TVARSTIS SU CINKO GELIU / NOBAZINCULAST® / NOBAVARO®
- LV** **Lietošanas instrukcija – ar cinka gēlu piesūcināti pārsēji**
PĀRSĒJS ZINKGELBINDE / NOBAZINCULAST® / NOBAVARO®
- ET** **Kasutusalaane teave**
TSINKGEELSIDE / NOBAZINCULAST® / NOBAVARO®
- HU** **Használati útmutató**
CINK GÉL KÖTÉS / NOBAZINCULAST® / NOBAVARO®
- RO** **Informații de utilizare**
FEȘELOR CU GEL PE BAZĂ DE ZINC / NOBAZINCULAST® / NOBAVARO®

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Product description:

Zinc gel bandages are gelatine-free half-rigid bandages characterized by an excellent thermal stability. They may consist of rigid (ZINC GEL BANDAGE) or lengthwise elastic (NOBAZINCULAST®) or crosswise and lengthwise elastic (NOBAVARO®) fabric.

Functional principle:

Zinc gel bandages are applied in moist condition, then harden and form a half-rigid shell with a compressing, supporting and stabilizing effect. NOBAZINCULAST® is lengthwise elastic before it hardens, while NOBAVARO® is lengthwise and crosswise elastic. These bandages are therefore easier to apply, without the need of cutting. Owing to the share of water contained, the bandages have a pleasant cooling effect (evaporation cooling), which is most remarkable with NOBAVARO® as this bandage contains a higher share of water. NOBAVARO® is hence particularly suited as a trauma bandage. As the viscosity of the gel is largely unaffected by the temperature, the trauma bandage can be stored in a cool place, for example the refrigerator, to increase the cooling effect.

Types of application:

ZINKGELBINDEN (rigid fabric)

NOBAZINCULAST® (fabric with longitudinal elasticity)

NOBAVARO (fabric with longitudinal and lateral elasticity)

Composition:

Zinc gel consisting of water, glycerin, zinc oxide, calcium alginate, sodium carboxymethyl cellulose, polyethylene glycol, 2-phenoxyethanol, methylparaben. Fabric material made of cotton (ZINKGELBINDE) resp. viscose and polyamide (NOBAZINCULAST® and NOBAVARO®)

Intended use and areas of application:

Zinkgel bandages are used as compression bandages, or else as support and relief bandages. They are applied by medical professionals.

- As compression bandage in the case of phlebological diseases, such as chronic weakness of the veins including ulcus cruris venosum (not in the weeping stage), venous oedema, thrombophlebitis, varicophlebitis, phlebothrombosis, in the case of posttraumatic oedema, oedema after achillodynia, lymphoedema.
- As support bandage and relief bandage in the case of diseases or injuries of the supporting apparatus and brace, sports injuries such as sprains, luxations, overstretchings, contusions, tendovaginitis, irritations of tendons, continued treatment after P.O.P. dressings.

Application instructions:

- When using rigid bandages please take care that no wrinkles are winded, for example at the joints. To avoid wrinkles, cut off the bandage several times and reposition it. It is also recommended to notch the edges of the bandage. The rest of the bandage can be applied without any problems.
- The advantage of the elastic bandages (NOBAZINCULAST® and NOBAVARO®) is the simpler application technique as they adapt to the respective surface and as they can also be applied on the joints without problems.
- Zinc gel bandages are winded in a circular way. They may be used directly on the skin.
- Open wounds should be covered with a suitable sterile wound pad before applying the bandage.
- When applying the zinc gel bandage, a good stability is achieved by winding four layers. If covered wounds exist, two additional layers should be applied.
- To avoid the soiling of clothes etc. NOBAFROTT® elastic terry cloth tube may be used as cover after the application of the zinc gel bandage. Make sure the zinc gel bandage is completely covered.
- If the dressing adheres too strongly to the skin when removing it, it is recommended to use hand-warm water to loosen the zinc gel.

Contra-indications:

- Zinc gel bandages may not be used with peripheric, arterial occlusions, decompensated cardiac insufficiency, septic phlebitis, phlegmasia coerulea dolens and on superinfected wounds.
- Do not apply the product in the case of a known allergy to the material.

Warnings and Precautions:

- Zinc bandages do not replace a medically indicated additional thrombosis prophylaxis.
- The reuse of single use devices may lead to a microbiological hazard. Treatment for reuse may adversely affect the performance of the product.
- Should serious incidents occur in relation to the device, they are, according to MDR (EU) 2017/745, to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Storage and transport:

Store in a horizontal position, do not store or transport at more than 35°C

Waste disposal:

In accordance with the requirements of the applicable local legislation and with the standards for the prevention of infections

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