

Specifications
Romed Elderly Care
Bath chairs
REF: BAT-01



Medical device conform MDR (EU) 2017/745

Introduction of the product

Romed Elderly Care is a product series with products specially designed for elderly people. These products give you a great comfort and will increase your independency.

Specifications

Romed bath chair

#	BAT-01
Specifications	PE back and seat with dripping holes, anti-skidding pad, easy to take and use
Material	Powder coated aluminum frame
Height	Adjustable
Tubing dimensions	Diameter=25.4 mm, Thickness=1.2 mm
Maximum user weight	110 kg
Net weight / gross weight	± 2 kg / ± 3

Specifications
Romed Elderly Care
Bath chair without back rest
REF: BAT-02



Medical device conform MDR (EU) 2017/745

Introduction of the product

Romed Elderly Care is a product series with products specially designed for elderly people. These products give you a great comfort and will increase your independency.

Specifications

Romed bath chair

#	BAT-02
Specifications	PE seat with dripping holes, anti-skidding pad, easy to take and use
Material	Powder coated aluminum frame
Height	Adjustable
Tubing dimensions	Diameter=25.4 mm, Thickness=1.2 mm
Maximum user weight	110 kg
Net weight / gross weight	± 2 kg / ± 3