

SPECIFICATIONS
ROMED ELECTRICAL WHEELCHAIRS
REF: WHE-POWER



Applicable harmonized standard: EN 12184:2009
Medical device conform MDR (EU) 2017/745

MATERIALS

Frame	Aluminium
Seat back cushion	Nylon

SIZES

Length	1040mm
Width	620mm
Height	900mm
Seat Width	425mm
Seat Height	485mm
Seat Depth	435mm
Armrest Hight	210mm
Armrest distance	465mm
Back Height	405mm
Front Wheel	8" (200mm) (Pneumatic)
Rear Wheel	22" (540mm) (Pneumatic)
Front Wheel Distance	500mm (Folding Width 420mm)
Rear Wheel Distance	555mm

MOTOR SPECIFICATIONS

Buzz Sound Level	Indoor 65dB/Outdoors 75dB
Motor SPEC	24V130W 3000RPM Tooth Proportion 49.41:1
Battery SPEC	24V18AH Lithium Battery
Control	Domestic Made
Charger	29.4V2A Charger (Lithium)
Endurance	More than 20Km
Climbing	Dynamic 6°, Static 9°

PROPERTIES

Fixed Armrest, Fix footrest, Foldable back, Aluminium Frame, Rear Wheel Motor Drive

Weight capacity: 100 Kg

Net weight: 29 Kg

PACKING

Per piece packed in an export carton