

**SPECIFICATIONS**  
**ROMED ELDERLY CARE**  
**MANUAL WHEELCHAIR DYNAMIC**  
**REF: WHE-02- BLACK**



**Reference number**  
WHE-02-BLACK

**Colour**  
Black

General specifications

Powder coated aluminum frame, Fixed full length armrest, PVC armpad, Swing away footrest with leg support belt, Backrest half folding, Hand strap, Safety belt, Nylon upholstery, Aluminum front fork, Plastic footplate

Seat width

48 cm

Wheels

8" solid PVC rear wheel and 8" solid PVC front wheel

Tubing material

High strength aluminum

Tubing diameter  
Tubing thickness

22.2 mm (25.4 mm for cross bar)  
2.0 mm (2.0 mm for cross bar)

Maximum user weight

110 kg

Net weight

± 11 kg

**Packing**

Per piece in an export carton