

SPECIFICATIONS

ROMED ANEROID SPHYGMOMANOMETERS

REF: BPM-AN-50



Applicable harmonized standard: EN ISO 81060-1:2012
Medical device conform MDD 93/42/EEC

Aneroid sphygmomanometer with adult sized cuff, with separate meter and bulb. A diagnostic apparatus for taking human blood pressure. Small in size and easy to carry. The diagnostic efficiency of the apparatus is similar to that of the mercurial type.

Specifications

Measuring range : 0-300 mmHg
Accuracy : ± 2.0 mmHg

Cuff size : 50 x 14 cm, one-size-fits-all for adults
Cuff material : Nylon with internal latex bladder and Velcro closure system, washable
Cuff colour : Blue, with artery indicator and gauge holder

Inflation bulb : Manual, with manual air release valve
Tubes length : 50 cm
Bulb, tube material : Rubber, black

Meter : Gauge to allow recalibration

Precautions

- When using this apparatus, do not let the pointer exceed 300 mmHg by inflation, because if the internal spring exceeds the elastic limit, accuracy can be affected.
- Avoid violent shocks in order to keep the reading accurate.
- Put the meter in the black vinyl carrying case after use and always keep it clean.

Packing

Per piece in an inner box, 50 pieces in an export carton.