

**SPECIFICATIONS OF
ROMED BLOOD LANCETS – STERILE DISPOSABLE
REF: BL-100TB**



Medical device conform MDD 93/42/EEC

Colour of plastic cover	: Blue
Cover material	: Polyethylene
Needle material	: Stainless steel ICr ₁₈ Ni ₉
Length of blood lancet (incl. cap on stitch point)	: 33mm
Thickness of stitch point	: 0.15 – 0.16 mm
Height of wrapper	: 6.2 ± 0.2 mm
Length of wrapper with lancet	: 23 mm
Lancet	: 3.2 ± 0.3 mm
Lancet needle size	: 23G
Sterilisation method	: Gamma Irradiation
Expiry date:	: 5 years after sterilisation
Dose	: 19.0 kGy; 39.0 kGy

Can be used with blood lancet pen (BL-PEN)

Packing

100 pieces in a dispenser box, 50 dispenser boxes in an inner box, 4 inner boxes in an export carton