

SPECIFICATIONS

ROMED BLOOD LANCETS - STERILE

REF: BL-200



Medical device conform MDD 93/42/EEC

Material of product	:	Stainless steel X ₆ Cr ₁₇
Length of blood lancet (incl. of stitch point)	:	39 ± 0.5 mm
Shape of lancet	:	Tribevel
Sharpness	:	≤0.8N
Lancet needle size	:	23G
Thickness of stitch point	:	0.15 – 0.16 mm
Sterilisation method	:	Gamma irradiation
Expiry date	:	5 years after sterilization

Packing: Per piece sterile in a paper pouch, 200 pcs in an inner box, 200 inner boxes (= total 40.000 pcs) in an export carton.