

## SPECIFICATIONS

### ROMED POLYETHYLENE GLOVES

REF: 65/65-LW

Reference number:	65/65-LW
Material:	100% high density polyethylene (HDPE), embossed
Colour:	Transparent
Sterility:	Non-sterile

#### DIMENSIONS

Contact thickness	0.011 mm
Length	30 cm ( $\pm 5\%$ )
Width	26 cm ( $\pm 5\%$ )
Weight	0.8 gram ( $\pm 5\%$ )

Packing: 100 pieces in a polybag, 100 polybags (=10.000 pieces) in an export carton