

SPECIFICATIONS
ROMED BLOOD LANCET PEN
REF: BL-PEN



Medical device conform MDR (EU) 2017/745

This product is used to puncture fingertip for blood sampling in hospitals and by diabetics.

Basic size	:	125 mm ± 1.75 mm
Material of product	:	Acrylonitrile butadiene styrene (ABS) and Polyoxymethylene (POM)
Elasticity of spring	:	65 Mn

Easy to open en close the cover, can be fit with lancet in right position.

Easy to adjust the dial for depth setting from 1 to 5.

At position 1, can inject through 4-6 ply foil, each with thickness 0.1 mm

At position 5, it can be injected through 16-18 ply of same foil.

Easy and safe on injection area, reliable to control and automatic withdrawal.

Easy and reliable to dispose the used lancet.

Can be used with Romed blood lancets BL-100TB.

Packing: Per piece in a box, 100 boxes in an export carton.