

SPECIFICATIONS

ROMED STETHOSCOPES

REF: STETH-SPL, STETH-DBL



Medical device conform MDD 93/42/EEC

Stethoscopes double head

Item no: STETH-DBL (double head)

- Anodized aluminum (alloy of zinc, aluminum and copper) rotating chest piece with built-in black rubber non-chill ring
- One-piece 'Y' tubing provides good transmission
- Length tube 53 cm
- Diaphragm diameter 40 mm, Bell diameter 20 mm
- External binaural spring
- One pair extra black rubber mushroom ear tips included for extra inquiry

Stethoscopes flat chest piece / single head

Item no: STETH-SPL (flat chest piece)

- Anodized aluminum (alloy of zinc, aluminum and copper) chest piece
- One-piece 'Y' tubing provides good transmission
- Length tube 53 cm
- Diaphragm diameter 40 mm
- External binaural spring
- One pair extra black rubber mushroom ear tips included for extra enquiry

Packing

Per piece in an inner box, 100 pieces in an export carton