

SPECIFICATIONS
ROMED LEG BAGS 750ML
WITH NON RETURN VALVE, 15CM TUBE, STERILE
REF: LB750MLNRV



Applicable harmonized standard: EN ISO 8669-2:1996
Medical device conform MDD 93/42/EEC

Bag: PVC
Softener: Soybean oil
Details: Plain bag without printing details
Thickness: 0.10mm in front and 0.10mm in back
Valve: Non-return valve inclusive
Bottom outlet: Push-pull valve (PVC and PE)
Tube length: 15cm
OD tube: 6.5mm

Universal fitting for connection with catheters.

White PVC connector and blue PE cap

Urine leg bags contain phthalates (DEHP).

Sterilization: EtO gas

Validity: 5 years after date of sterilization

Packing: Per piece sterile in a bag, 10 bags in a polybag and 250 pieces in an export carton

VAN OOSTVEEN MEDICAL B.V. – ROMED HOLLAND, Herenweg 269, 3648 CH Wilnis, NL
www.romed.nl Tel.: +31 0297 282101, fax: +31 0297 288316, e-mail: info@romed.nl

2017-08