

SPECIFICATIONS
ROMED ELDERLY CARE
MANUAL WHEELCHAIR DYNAMIC
REF: WHE-01-RED/BLUE/BLACK



Applicable harmonized standard: EN 12183:2014
Medical device conform MDR (EU) 2017/745

Reference number

WHE-01-RED
WHE-01-BLUE
WHE-01-BLACK

Colour

Red
Blue
Black

General specifications

Powder coated steel frame, Flip-up armrest. Swing away footrest with leg support belt, nylon upholstery, steel handrim and plastic footplate

Seat width

46 cm

Wheels

24" solid rubber rear wheel and 8" front caster

Tubing material

High strength steel

Tubing diameter
Tubing thickness

22.2 mm (25.4 mm for cross bar)
1.2 mm (1.6 mm for cross bar)

Maximum user weight

110 kg

Net weight

± 17 kg

Packing

Per piece in an export carton