SPECIFICATIONS ROMED URINE BAGS 1500CC AND 2000CC, STANDARD REF: UB1.5LSTAND, UB2LSTAND

ROMED URINE BAGS 1500CC AND 2000CC WITH NON RETURN VALVE REF: UB1.5LNRV, UB2LNRV

Applicable harmonized standard: EN ISO 8669-2:1996 Medical device conform MDR (EU) 2017/745

Standard: UB1.5LSTAND / UB2LSTAND

Urine collection bag 1500cc/2000cc manufactured in translucent PVC with scaled measurements from 100cc to 1500cc/2000cc. Tube length 90 cm. Translucent tube hermetically sealed to the bag and fitted with a

conical connector at the end adaptable to the patient's drain.

With non-return valve: UB1.5LNRV / UB2LNRV

Urine collection bag 1500cc/2000cc manufactured in translucent PVC with scaled measurements from 100cc to 1500cc/2000cc. Tube length 90 cm.

Incorporates a double plate in the upper part acting as a valve to avoid the return of fluids and gases. Translucent tube hermetically sealed to the bag and fitted with a conical connector at the end adaptable to the patient's drain.

Component	Parameter	Specification
Bag	Foil Soldering Capacity Softener	PVC High frequency 1500cc or 2000cc Soybean oil
Non return valve	Material	PVC
Tube	Length Diameter Material	90 cm 0.5 cm PVC, flexible
Connector	Material	PVC, white
Тар	Material	LDPE

TECHNICAL SPECIFICATIONS

Urine bags contain phthalates (DEHP) Universal fitting for connection with catheters

Packing

UB1.5LSTAND / UB2LSTAND: 10 pieces in a bag, 250 pieces in an export carton UB1.5LNRV / UB2LNRV: 10 pieces in a bag, 250 pieces in an export carton

VAN OOSTVEEN MEDICAL B.V. – ROMED HOLLAND, Herenweg 269, 3648 CH Wilnis, NL www.romed.nl Tel.: +31 0297 282101, fax: +31 0297 288316, e-mail: info@romed.nl