

# SPECIFICATIONS

## ROMED DIGITAL BABY SCALE

### REF: BS-003



BS-003 Romed baby scale is a high precision digital weighing scale.

Maximum weight : 20.000 kg  
Precision : 0.005 g

Accuracy range 0-5kg : 0.025 kg  
Accuracy range 5-10kg : 0.030 kg  
Accuracy range 10-20kg : 0.040 kg

Units : kg and lbs

#### CHARACTERISTICS

LCD display  
Hold function for recording weight of moving child  
Tare function  
Automatic zero and switch off  
Low battery and overload indication

Power : 2x AAA batteries (included)  
Material : ABS plastic  
Dimension scale : 56.0 x 33.0 x 4.6 cm  
Net weight scale : 1.68 kg

#### LEGISLATION

Conform Electromagnetic Compatibility (EMC) Directive 2014/30/EU  
Conform RoHS Directive 2011/65/EU

#### PACKING

Per piece in a box and 4 pieces in an export carton

VAN OOSTVEEN MEDICAL B.V. – ROMED HOLLAND, Herenweg 269, 3648 CH Wilnis, NL  
[www.romed.nl](http://www.romed.nl) Tel.: +31 0297 282101, fax: +31 0297 288316, e-mail: [info@romed.nl](mailto:info@romed.nl)

2018-01