

**SPECIFICATIONS**  
**ROMED PENLIGHTS DELUXE**  
**REF: PEN-DLB**



*Medical device conform MDR (EU) 2017/745*

**PROPERTIES**

Material : Aluminium  
Colour : Blue  
Light : LED  
Battery : 2x AAA, batteries can be replaced

With pinpoint, high intensity white light bulb, 3 Volts

For medical use.

**DIMENSIONS**

Length : 12.7 cm  
Diameter : 1.3 cm  
Weight : 47 g

**PACKING**

Packed per piece in a suede etui, 25 pcs in a dispenser box, 4 dispenser boxes (=100 pieces) in an export carton