

SPECIFICATIONS
ROMED DISPOSABLE HOSPITAL RAZORS
REF: RAZORS



Medical device conform MDR (EU) 2017/745

DESCRIPTION

Type: Single edged one-blade razor with plastic handle, non sterile

MATERIALS

Handle: Polyethylene
 Blades: Stainless steel (6Cr13) with single edge coated
 3-ply coating: Cr + Ceramal + PTFE Teflon

CHARACTERISTICS

Blades

Length	36.8 ± 0.05mm
Width	7.14 ± 0.01mm
Thickness	0.1 mm
Sharpness	≤ 20N
Hardness	590-650 HV

Handle

Length	98.5 ± 0.1mm
--------	--------------

PACKING

100 pieces in a dispenser box, 12 boxes (=1200 pieces) in an export carton