

SPECIFICATIONS

ROMED DIGITAL THERMOMETER WATERPROOF

REF: THERM-DIG-WP



Applicable harmonized standard: EN ISO 80601-2-56:2012
Medical device conform MDD 93/42/EEC

Measurement range:	32.0~42.9°C (Display Lo°C when the temperature is below 32.0°C Display Hi°C when the temperature is above 42.9°C)
Accuracy:	±0.1°C (35.5~42°C) ±0.2°C (<35.5°C or >42°C)
Display mode:	Liquid crystal display, 3.5 digits
Power consumption:	0.15mW under operating state
Battery:	1.5V button battery (LR41), replaceable by either an LR41 or SR41 battery
Battery life (in use):	About 200 hours under continuous operation
Battery life (in storage):	The battery will maintain ≥90% of the capacity after one year storage at 20±2°C and 55±20% RH
Thermometer shelf life:	The thermometer does not have a defined shelf life
Dimensions:	128 x 18 x 11 mm
Net Weight:	About 10 grams
Alarm:	About 10 seconds sound after operation

Scale: Celsius

Dispose batteries as chemical waste

Packing

Packed per piece in a plastic case, then in a small box, 20 pieces in an inner box, 240 pieces in an export carton