

SPECIFICATIONS

ROMED DIGITAL THERMOMETER

REF: THERM-DIG + THERM-FLEX



Applicable harmonized standard: EN ISO 80601-2-56:2012
Medical device conform MDD 93/42/EEC

THERM-DIG: Normal type
 THERM-FLEX: Flexible tip type

Measurement range: 32~42°C
 (Display L°C when the temperature is below 32°C
 Display H°C when the temperature is above 42°C)

Accuracy: ±0.1°C (35~40°C)
 ±0.2°C (<35°C or >42°C)

Display mode: Liquid crystal display, 3.5 digits

Power consumption: 0.15mW under operating state

Battery: 1.5V button battery (LR-41)

Battery life: About 100 hours under continuous operation

Dimensions: 129 x 22 x 14mm (including transparent box)

Net Weight: ± 11 grams

Alarm: ± 4 seconds sound after operation

Scale: Celsius

Dispose batteries as chemical waste

Packing: packed per piece in a plastic case, then in a small box, 20 pieces in an inner box, 200 pieces in an export carton