

SPECIFICATIONS

ROMED LATEX FINGER COTS

REF: FC-S, FC-M, FC-L, FC-XL



Medical device conform MDR (EU) 2017/745

Material	:	Natural Rubber Latex DRC 60%
		Pre powdered
Average weight per piece:	Small	: 0.42 ± 0.02 g
	Medium	: 0.47 ± 0.02 g
	Large	: 0.52 ± 0.02 g
	X-large	: 0.59 ± 0.02 g
Diameter of Bead	:	1.20 ± 0.20 mm
Average length	:	70 ± 3 mm
Tensile strength at Break (ASTM D 412)	:	≥ 15 MPa
Elongation at break	:	$\geq 700\%$
Thickness	:	0.09 ± 0.03 mm
Resistance level for Anti-Static types	:	$10^9 - 10^{11} \Omega/\text{cm}^3$
Thickness at Tip	:	≥ 0.11 mm
Width when laid flat	Small	: 23 ± 2 mm
	Medium	: 27 ± 2 mm
	Large	: 29 ± 2 mm
	X-large	: 33 ± 2 mm

Caution: This product contains natural rubber latex which may cause allergic reactions

Packing: 100 pieces in a plastic bag, 250 bags (=25.000 pieces) in an export carton

VAN OOSTVEEN MEDICAL B.V. – ROMED HOLLAND, Herenweg 269, 3648 CH Wilnis, NL
www.romed.nl Tel.: +31 0297 282101, fax: +31 0297 288316, e-mail: info@romed.nl

2021-05