

SPECIFICATIONS

ROMED 3-WAY STOP COCK

REF: 3WSC



Applicable harmonized standard: EN 20594-1:1993
 Medical device conform MDD 93/42/EEC

3-Way Stop Cocks are designed to provide an additional infusion opportunity between in line and venepuncture. Fully transparent polycarbonate body for visualisation of flow.
 6% Luer taper for compatibility with any international product. Arrow indication marks for indicating direction of flow.

Component	Raw material
Channel	Polycarbonate
Handle	Polyethylene, HDPE 50 MA 180
Screw cap	Polyvinyl chloride (PVC)
Luer lock	Polyethylene
Cock cap	Polypropylene
Packaging	Medical grade paper/Tyvek with standard plastic film

This product does not contain phthalates.

Characteristics

- Rotating dead space for accurate drug administration
- Continuous flow channels
- Tight covers for all ends
- Rotation allowed for 360°
- Pressure tested pneumatically and hydrostatically for up to 60 psi

Sterilization

Sterilized by ethylene oxide

Packing

Sterile per piece, 50 pieces in an inner box and 20 x 50 = 1000 pieces in an export carton