

**SPECIFICATIONS**  
**ROMED SHOE COVERS**  
**REF:SC-2000**



*Medical device conform MDR (EU) 2017/745*

Material : Cast polyethylene film (embossed)

Colour : Blue

**Dimensions**

Contact thickness : 0.050 mm

Length : 360 mm ( $\pm 5\%$ )

Width : 150 mm ( $\pm 5\%$ )

Weight : 3.2 gram per piece

Shoe covers are latex free.

These products are medical grade personal protection disposable devices and therefore are also suitable for use in the food industry.

**Packing**

10 pieces per bundle, 10 bundles in a clear polybag, 2000 pieces in an export carton

VAN OOSTVEEN MEDICAL B.V. – ROMED HOLLAND, Herenweg 269, 3648 CH Wilnis, NL  
[www.romed.nl](http://www.romed.nl) Tel.: +31 0297 282101, fax: +31 0297 288316, e-mail: [info@romed.nl](mailto:info@romed.nl)

2021-05