

SPECIFICATIONS
ROMED WEEKLY MEDICATION DISPENSER
REF: PD-WEEK



Medical device conform MDR (EU) 2017/745

The Romed Weekly Medication Dispenser is a tool for the daily distribution of medicines. The product is a holder with seven removable parts for every day of the week. The individual parts have a clear sliding lid with four compartments for the saving of pills per part of the day.

Different days of the week and day parts are indicated in English and German.

DIMENSIONS

Length : 115 mm
Width : 43 mm
Height : 127 mm

Length day part : 107 mm
Width day part : 34 mm
Height day part : 16 mm

Weight : 171 gram

PACKING

Packed per piece, 60 pieces in an export carton