

ROMED COTTON TIPPED APPLICATOR, STERILE, 150MM

SPECIFICATIONS

REF: ECO-CTA15-ST



Medical device conform MDR (EU) 2017/745

INTENDED USE

A wooden stick with a cotton head used to clean the surface of the skin. The highly absorbent cotton head has strong water absorption, and after absorbing the disinfectant, the disinfectant can wipe the skin evenly to achieve disinfection effect. This is a single use device.

MATERIALS

Stick: wood
Head: cotton

SIZES

150 mm

SPECIFICATIONS

Wooden stick, natural colour
Cotton head, white
Clean and smooth appearance, no burr.
Sterile

CAUTION AND WARNINGS


This is a single use device. Do not reuse.
Keep away from sunlight.
Please dispose of the product after use in compliance with local regulation.

SHELF LIFE:

5 years from the production date

PACKING

Per piece in a pouch, 100 pieces in a box, 50 boxes = 5000 pieces in an export carton.



GB: <https://www.romed.nl/en/products/hygiene-protection/cotton-tipped-swabs/cotton-tipped-applicators/groups/q+c+a+view>
NL: <https://www.romed.nl/nl/producten/hygiene-bescherming/mondverzorging/wattenstaafjes/groups/q+c+a+view>
DE: <https://www.romed.nl/de/produkte/hygiene-schutz/mundpflege/wattestabchen/groups/q+c+a+view>
FR: <https://www.romed.nl/fr/produits/hygiene-protection/soins-bucco-dentaires/applicateurs-en-coton/groups/q+c+a+view>
ES: <https://www.romed.nl/es/productos/higiene-y-proteccion/bastoncillos-para-el-cuidado-bucal/aplicadores-de-algodon/groups/q+c+a+view>
PT: <https://www.romed.nl/pt/productos/higiene-e-protecao/esfregacos-de-cuidado-oral/aplicadores-de-algodoao/groups/q+c+a+view>
IT: <https://www.romed.nl/it/prodotti/igiene-e-protezione/bastoncini-tampone-per-igiene-orale/applicatore-in-cotone/groups/q+c+a+view>

