

ROMED COTTON TIPPED APPLICATOR, NON-STERILE

SPECIFICATIONS

REF: ECO-CTA15, ECO-CTA7.5



Medical device conform MDR (EU) 2017/745

INTENDED USE

A wooden stick with a cotton head used to clean the surface of the skin. The highly absorbent cotton head has strong water absorption, and after absorbing the disinfectant, the disinfectant can wipe the skin evenly to achieve disinfection effect. This is a single use device.

MATERIALS

Stick:

ECO-CTA15	Wood, natural colour
ECO-CTA7.5	Paper, white

Head: cotton

SIZES

ECO-CTA15	150 mm
ECO-CTA7.5	75 mm

SPECIFICATIONS

ECO-CTA15	Single tipped
ECO-CTA7.5	Double tipped

Paper stick, white

Cotton head, white

Clean and smooth appearance, no burr.

Non-sterile

CAUTION AND WARNINGS

This is a single use device. Do not reuse.

Keep away from sunlight.

Please dispose of the product after use in compliance with local regulation.

SHELF LIFE:

3 years from the production date

PACKING

ECO-CTA15	100 pieces in a bag, 10 bags = 1000 pieces in a box, 10 boxes = 10000 pieces in an export carton.
ECO-CTA7.5	200 pieces in a box, 50 boxes = 10000 pieces in an export carton.



GB: <https://www.romed.nl/en/products/hygiene-protection/cotton-tipped-swabs/cotton-tipped-applicators/groups/q+c+a+view>
 NL: <https://www.romed.nl/nl/producten/hygiene-bescherming/mondverzorging/wattenstaafjes/groups/q+c+a+view>
 DE: <https://www.romed.nl/de/produkte/hygiene-schutz/mundpflege/wattestabchen/groups/q+c+a+view>
 FR: <https://www.romed.nl/fr/produits/hygiene-protection/soins-bucco-dentaires/applicateurs-en-coton/groups/q+c+a+view>
 ES: <https://www.romed.nl/es/productos/higiene-y-proteccion/bastoncillos-para-el-cuidado-bucal/aplicadores-de-algodon/groups/q+c+a+view>
 PT: <https://www.romed.nl/pt/produtos/higiene-e-protecao/esfregacos-de-cuidado-oral/aplicadores-de-algodoao/groups/q+c+a+view>
 IT: <https://www.romed.nl/it/prodotti/igiene-e-protezione/bastoncini-tampone-per-igiene-orale/applicatore-in-cotone/groups/q+c+a+view>