


ROMED IV DRESSING, waterproof, non-woven

REF: IV1200

SPECIFICATIONS



Medical device conform MDR (EU) 2017/745

INTENDED USE

The Romed IV Dressing is intended to be used for securing and covering IV catheters.

MATERIALS

Layer 1, cover paper, 55g/m²

Layer 2, viscose fibre, 160g/m² + isolation film

Layer 3, non-woven fiber, 50g/m², with solvent pressure sensitive adhesive

SIZES

One size: 6 x 8 cm

SPECIFICATIONS

- Cannula Dressing
- Superior absorption
- Simple to apply
- Gentle and moldable

INSTRUCTIONS FOR USE

- Apply to clean, dry skin
- Discontinue use if persistent skin irritation exists
- Store in a cool dry place
- Sterile unless package is opened

SHELF LIFE: 5 years from the production date

PACKING

Per piece is sterile packaging, 50 pieces per box, 1200 pieces per export carton

