

SPECIFICATIONS

ROMED SPECIAL CONDOMS (CLASSIC, COLOURED, EXTRA THIN, EXTRA STRONG, RIBBED & DOTTED, FLAVOURED)

REF: CON3P40, CON-COL, CON-THN, CON-STR, CON-TEX, CON-FLA



Applicable harmonized standard: EN ISO 4074:2015
Medical device conform MDD 93/42/EEC

Latex prophylactics from pure high strength, natural rubber latex in accordance with the international standards, and free from cosmetic defects. Inspection : (100%) electronically

Caution: This product contains natural rubber latex which may cause allergic reactions.

Condoms do not contain nitrosamines. They do contain traces of NDMA 90, but within the specifications as stipulated in directive 93/11/EEC.

Classic (CON3P40) (package of three pieces)

Type	:	Straight sided, smooth surface
Style	:	Reservoir end
Surface treatment	:	Lubricated with silicone oil
Amount of lubrication	:	≤ 400 mg
Dressing material	:	Corn starch powder, 40-55 mg
Colour	:	Pink
Length	:	180 ± 10 mm
Width	:	52 mm ± 2 mm
Thickness	:	0.065 ~ 0.08 mm

COLOURED (CON-COL) (package of three pieces)

Type	:	Straight sided, smooth surface
Style	:	Reservoir end
Surface treatment	:	Lubricated with silicone oil
Amount of lubrication	:	≤ 400 mg
Dressing material	:	Corn starch powder, 40-55 mg
Colour	:	Red, pink and yellow
Length	:	180 ± 10 mm
Width	:	52 mm ± 2 mm
Thickness	:	0.065 ~ 0.08 mm

FLAVOURED (CON-FLA) (package of three pieces)

Type	:	Straight sided, smooth surface
Style	:	Reservoir end
Surface treatment	:	Lubricated with silicone oil premixed with flavour
Amount of lubrication	:	≤ 400 mg
Dressing material	:	Corn starch powder, 40-55 mg
Colour	:	Red, orange and yellow
Length	:	180 ± 10 mm
Width	:	52 mm ± 2 mm
Thickness	:	0.065 ~ 0.08 mm

RIBBED & DOTTED (CON-TEX) (package of three pieces)

Type	:	Straight sided, one fully dotted, one fully ribbed, one fully dotted and ribbed
Style	:	Reservoir end
Surface treatment	:	Lubricated with silicone oil premixed with flavour
Amount of lubrication	:	≤ 400 mg
Dressing material	:	Corn starch powder, 40-55 mg
Colour	:	Natural
Length	:	180 ± 10 mm
Width	:	52 mm ± 2 mm
Thickness	:	0.065 ~ 0.08 mm

EXTRA THIN (CON-THN) (package of three pieces)

Type	:	Straight sided, smooth surface
Style	:	Reservoir end
Surface treatment	:	Lubricated with silicone oil
Amount of lubrication	:	≤ 400 mg
Dressing material	:	Corn starch powder, 40-55 mg
Colour	:	Natural
Length	:	180 ± 10 mm
Width	:	52 mm ± 2 mm
Thickness	:	0.060 ± 0.015 mm

EXTRA STRONG (CON-STR) (package of three pieces)

Type	:	Straight sided, smooth surface
Style	:	Reservoir end
Surface treatment	:	Lubricated with silicone oil
Amount of lubrication	:	≤ 400 mg
Dressing material	:	Corn starch powder, 40-55 mg
Colour	:	Natural
Length	:	180 ± 10 mm
Width	:	52 mm ± 2 mm
Thickness	:	> 0.100 mm

Physical properties

Classic, coloured, flavoured, ribbed & dotted, extra thin

Volume at burst: ≥ 18.0 dm³

Pressure at burst: ≥ 1.00 kPa

Extra strong condoms:

Elongation at break: ≥ 700%

Force at break: ≥ 100 N

Volume at burst: ≥ 18.0 dm³

Pressure at burst: ≥ 2.00 kPa

Packing

A strip of 3 pieces in a small box, 12 small boxes in a display box, 4 display boxes in an inner box, 30 inner boxes in an export carton