

# SPECIFICATIONS

## ROMED NON-WOVEN FACE MASK, 3 PLY, WITH EARLOOPS

### REF: MASK-L



Medical device conform MDR (EU) 2017/745

<b>Applicable standard</b>	EN 14683:2019, Type IIR – mentioned on the packaging
<b>Appearance</b>	Non-Woven Face Mask, 3PLY 2 elastic loops, pleats and nose clip
<b>First Ply</b>	Non-Woven Cloth, Spunbond PP, White, 25 ±2 g/m <sup>2</sup>
<b>Second Ply</b>	Filter Paper, PP, white, 22 ±2 g/m <sup>2</sup>
<b>Third Ply</b>	Non-Woven Cloth, Spunbond PP, blue, 18 ±2 g/m <sup>2</sup>
<b>Size</b>	Length: 18 ±0.5cm Width unfolded: 16.5 ±0.5cm Width folded: 9.5 ±0.5cm
<b>Bacterial Filtration Efficiency (BFE)</b>	≥99%
<b>Delta P (Mean Particle Size 3.0 µm)</b>	<49 Pa/m <sup>2</sup>
<b>Fluid resistance</b>	16 kPa (120 mmHg)

#### General properties:

- Latex free
- Biocompatible
- Flammability Type A, Class 1
- Lint free
- Hypoallergenic
- Breathable
- Tear resistant
- Fluid resistant layer
- Integrated flexible nose clip (nickel free)
- Ultrasonic bonding in the areas where the mask is connected with the ties or loops in order to ensure impermeability to liquids.

**Packing:** 50 pieces in a box, 20 boxes (= 1000 pieces) in an export carton