

SPECIFICATIONS
ROMED ELASTIC BANDAGES
REF: EB-



Medical device conform MDR (EU) 2017/745

Materials

57% Cotton
39% PBT (polybutylene terephthalate)
4% Elastane

Clip: Aluminium + elastic ribbon with 30% latex content, no skin contact

Product is PVC and DEHP-free

Material properties

Weight: 68 g/m²
Elasticity: ≥ 160%

Sizes

Article number	Size
EB8x5	8 cm x 5 m
EB10x5	10 cm x 5 m
EB12x5	12 cm x 5 m
EB15x5	15 cm x 5 m

Packaging

EB8x5 Packed per roll in an inner box, 250 rolls in export carton
EB10x5 Packed per roll in an inner box, 200 rolls in export carton
EB12x5 Packed per roll in an inner box, 200 rolls in export carton
EB15x5 Packed per roll in an inner box, 100 rolls in export carton