

**SPECIFICATIONS**  
**ROMED SILK TAPE**  
**REF: ST-1.25, ST-2.5, ST-5**



*Medical device conform MDR (EU) 2017/745*

**MATERIALS**

Silk tape is made of silk fabric, coated with hot melt glue on the base fabric

Silk tape is also called micropore tape.

**INTENDED USE**

Silk tape is commonly used to secure electrodes/angiocath/heparin block/tubings/dressing materials to the patient's skin.

**SIZES**

Reference number	Width	Length
ST-1.25	1.25 ±0.1 cm	5.0 m
ST-2.5	2.5 ±0.1 cm	5.0 m
ST-5	5.0 ±0.1 cm	5.0 m

**CAUTION AND WARNINGS**

- This is a single use product, it should not be reused
- This product should not be used on irritated skin
- This product should not be placed directly on a wound

**Shelf life:** 5 years from the production date

**PACKING**

Reference number	Box	Carton
ST-1.25	24 rolls	720 rolls
ST-2.5	12 rolls	360 rolls
ST-5	6 rolls	180 rolls