SPECIFICATIONS ROMED ALCOHOL PREPS REF: PREP-2000

IDENTIFICATION

Product name: ROMED Alcohol preps Ref: PREP-2000

PURPOSE

For cleaning purpose.

COMPOSITION

Prep: non-woven material (spunlace, mixture of polyester and viscose) Weight: 50 g/m² Thickness: 0.35 mm, 2-ply, Dimensions prep: 65 x 30 mm \pm 2 mm Active ingredient: 70% isopropyl alcohol Inactive ingredient: water Liquid 0.4 – 0.6 ml per package Alcohol contents \pm 0.02 g/cm², which is \pm 0.45 per piece

Foil package: paper, aluminium foil and plastic. weight single pack: 73 g/m², thickness: 0.06 mm

PRECAUTIONS

- Avoid contact with eyes. If contacted, flush eyes with copious quantities of water.
- If taken in, DO NOT induce vomiting, but obtain medical attention immediately.
- Stop use if irritation and redness develop on the skin. Ask a doctor if condition persists for more than 72 hours

SPECIAL NOTES

Single use, discard after use. Keep at room temperature. For external use only. Flammable, keep away from fire or flame. Keep out of reach of children.

Packing: 100 pieces, each packed in paper pouch, in an inner box, 50 inner boxes in an export carton.