

SPECIFICATIONS
ROMED POLYETHYLENE COAT APRONS
REF: 21/60-CA



Medical device conform MDR (EU) 2017/745

REFERENCE

21/60-CA

Blue, packed per piece in a polybag

MATERIAL

100% Chlorinated polyethylene (CPE) film, embossed

DIMENSIONS

Length	116 cm (\pm 5%)
Waist length	98 cm (\pm 5%)
Length Arm-Arm	170 cm (\pm 5%)
Weight	36.5 g (\pm 5%)

PACKING

Packed per piece in a bag, 10 bags in an inner box, 30 inner boxes (=300 pieces) in an export carton.