

SPECIFICATIONS

ROMED STANDARD CONDOMS

REF: CON-500R



Applicable harmonized standard: EN ISO 4074:2015
Medical device conform MDD 93/42/EEC

Latex prophylactics from pure high strength, natural rubber latex in accordance with the international standards, and free from cosmetic defects. Inspection : (100%) electronically

Type : Straight sided, smooth surface
Style : Reservoir end
Surface treatment : Lubricated with silicone oil
Amount of lubrication : ≤ 400 mg
Colour : Pink

Dimensions :
* Length : 180 ± 10 mm
* Width : 52 ± 2 mm (measured at a distance of 80 ± 5 mm from the rim when laid flat)
* Wall thickness : $0.065 \sim 0.08$ mm
* Mass : $1.40 \sim 1.60$ g

Physical properties

Bursting volume : $\geq 18.0 \text{ dm}^3$
Pressure at burst : $\geq 1.00 \text{ kPa}$

Inspection : 100% electronically

Caution: This product contains natural rubber latex which may cause allergic reactions.

Condoms do not contain nitrosamines. They do contain traces of NDMA 90, but within the specifications as stipulated in directive 93/11/EEC.

Packing

CON500R: 48 strips of 3 pieces (144 pieces) in an inner box, 50 inner boxes (7.200 pieces) in an export carton.