

SPECIFICATIONS
ROMED GAUZE COMPRESSES NON-STERILE
REF: GCNS-050, GCNS-075, GCNS-100



Applicable harmonized standard: EN 14079:2003
Medical device conform MDR (EU) 2017/745

MATERIAL

Material: 100% Cotton
Layers: 8 Ply
Threads: 13

Non-Sterile

SIZES

Reference number	Width	Length
GCNS-050	5.0 cm	5.0 cm
GCNS-075	7.5 cm	7.5 cm
GCNS-100	10.0 cm	10.0 cm

PROPERTIES

Colour	White
Breaking load (Warp direction)	35N/500mm
Breaking load (Weft direction)	20N/500mm
Weight	16.5 g/m ²

PACKING

GCNS-050: 100 pieces in paper pack, 200 paper packages in an export carton
GCNS-075: 100 pieces in paper pack, 100 paper packages in an export carton
GCNS-100: 100 pieces in paper pack, 100 paper packages in an export carton