



# TOP BRUSH

## DISPOSABLE TOOTHBRUSHES

### CONFIGURATION OVERVIEW

**Type:**

- non-sterile
- individually wrapped
- pre-pasted

### PRODUCT VARIATIONS

Color Code	Y	O	P	L	FG	MB	B	WR	F	DG
Color	Yellow	Orange	Pink	Lilac	Fresh green	Med blue	Blue	Burgundy	Fuchsia	Green
										

**Art-No.**

ZABU2-AZ-Color Code

ZABU2-AZ-MIX (25 Y, 25 O, 25 DG, 25 B)  
 ZABU2-AZ-MIX1 (25 MB, 25 WR, 25 FG, 25 P)

## TECHNICAL SPECIFICATIONS

### Body:

Material	PP
Length	164 mm $\pm$ 1.0 mm
Width (head)	10.8 mm $\pm$ 0.3 mm
No. of tuft holes	39
Depth of tuft holes	3.2 mm $\pm$ 0.2 mm
Diameter of tuft holes	1.6 mm $\pm$ 0.05 mm

### Filament:

Material	PP
Type (rounding)	Flat
Diameter	0.18 mm $\pm$ 0.01 mm
Length (at highest point)	10 mm $\pm$ 0.5 mm
Stiffness	Medium
No. of filament/tuft	22 $\pm$ 2

### Staple:

Material	Aluminum wire
Length	2.00 mm $\pm$ 0.15 mm
Width	1.30 mm $\pm$ 0.10 mm
Thickness	0.25 mm $\pm$ 0.05 mm

### Paste solution:

Ingredients	Deionized water, alcohol, saccharin, xylitol, sodium, lauryl sulfate, menthol, peppermint oil
Color	Transparent

## PACKAGING

<b>Pieces/box:</b>	100
<b>Carton size:</b>	57x36.5x19.5 cm
<b>Carton/layers:</b>	4
<b>Carton/pallet:</b>	40
<b>Pallet size:</b>	80x120x210 cm

<b>Boxes/carton:</b>	10
<b>Carton weight:</b>	7.5 kg
<b>Layers/pallet:</b>	10
<b>Pallet weight:</b>	320 kg